ROBERT KOCH INSTITUT

KL FÜR POCKENVIREN

# Präanalytikhandbuch

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Diaz-Cánova D, Moens UL, **Brinkmann A, Nitsche A, Okeke MI (2022).** Genomic Sequencing and Analysis of a Novel Human Cowpox Virus With Mosaic Sequences From North America and Old World Orthopoxvirus. *Front Microbiol* 13: 868887.

Bourquain D, **Schrick L, Tischer BK, Osterrieder K, Schaade L, Nitsche A (2021).** Replication of cowpox virus in macrophages is dependent on the host range factor p28/N1R. *Virology J* 18 (1): 173.

Mavian C, López-Bueno A, Martín R, **Nitsche A, Alcamí A (2021).** Comparative pathogenesis, genomics and phylogeography of mousepox. *Viruses* 13 (6): 1146.

Wendt R, Tittelbach J, **Schrick L, Kellner N, Kalbitz S, Ruehe B, Michel J, Schliemann S, Elsner P, Lübbert C, Nitsche A (2021).** Generalized cowpox virus infection in an immunosuppressed patient. *Int J Infect Dis* 106 (May): 276–278.

Esparza J, **Nitsche A, Damaso CR (2020).** Investigations on the historical origin and evolution of the smallpox vaccine. *Gaceta Médica de Caracas* 128 (Supl 1): S88–S97.

**Brinkmann A**, Souza ARV, Esparza J, **Nitsche A**, Damaso CR (2020). Re-assembly of nineteenth-century smallpox vaccine genomes reveals the contemporaneous use of horsepox and horsepox-related viruses in the USA. *Genome Biol* 21 (1): 286.

Esparza J, Lederman S, **Nitsche A**, Damaso CR (2020). Early smallpox vaccine manufacturing in the United States: Introduction of the “animal vaccine” in 1870, establishment of “vaccine farms”, and the beginnings of the vaccine industry. *Vaccine* 28 (30): 4773–4779.

Patrono LV, Pléh K, Samuni L, Ulrich M, Röthemeier C, Sachse A, **Muschter S**, **Nitsche A**, Couacy-Hymann E, Boesch C, Wittig RM, Calvignac-Spencer S, Leendertz FH (2020). Monkeypox virus emergence in wild chimpanzees reveals distinct clinical outcomes and viral diversity. *Nat Microbiol* 5 (7): 955-965.

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Esparza J, **Nitsche A**, Damaso CR (2018). Beyond the myths: Novel findings for old paradigms in the history of the smallpox vaccine. *PLoS Pathog* 14(7):e1007082.

Esparza J, **Schrack L**, Damaso CR, **Nitsche A** (2017). Equination (inoculation of horsepox): an early alternative to vaccination (inoculation of cowpox) and the potential role of horsepox virus in the origin of the smallpox vaccine. *Vaccine* 35(52): 7222-7230.

Stagegaard J, Kurth A, Stern D, Dabrowski PW, Pocknell A, **Nitsche A**, **Schrack L** (2017). Seasonal recurrence of Cowpox virus outbreaks in captive cheetahs (*Acinonyx jubatus*). *PLoS One* 12(11): e0187089.

Grönemeyer LL, Baltzer A, Broekaert S, **Schrack L**, Möller L, **Nitsche A** et al. (2017). Generalised cowpox virus infection. *Lancet* 390(10104): 1769.

**Grossegeesse M**, Doellinger J, Haldemann B, Schaade L, **Nitsche A** (2017): A Next-Generation Sequencing Approach Uncovers Viral Transcripts Incorporated in Poxvirus Virions. *Viruses* 9 (10): pii: E296.

**Schrack L**, **Tausch SH**, Dabrowski PW, Damaso CR, Esparza J, **Nitsche A** (2017). An early American smallpox vaccine based on horsepox. *N Engl J Med* 377(15): 1491–1492.

Wibbelt G, **Tausch SH**, Dabrowski PW, Kershaw O, **Nitsche A**, **Schrack L** (2017). Berlin Squirrelpox Virus, a new poxvirus in red squirrels, Berlin, Germany. *Emerg Infect Dis* 23(10): 1726–1729.

Leendertz SAJ, Stern D, Theophil D, Anoh E, Mossoun A, Schubert G, Wiersma L, Akoua-Koffi C, Couacy-Hymann E, Muyembe-Tamfum JJ, Karhemere S, Pauly M, **Schrack L**, Leendertz FH, **Nitsche A** (2017). A Cross-sectional Serosurvey of anti-Orthopoxvirus Antibodies in Central and Western Africa. *Viruses* 9(10): pii: E278.

Stern D, Olson VA, Smith SK, Pietraszczyk M, **Miller L**, Miethe P, Dorner BG, **Nitsche A** (2017). Rapid and sensitive point-of-care detection of Orthopoxviruses by ABICAP immunofiltration. *Viol J* 13(1): 207.

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Stern D, Pauly D, Zydek M, Miller L, Piesker J, Laue M, Lisdat F, Dorner MB, Dorner BG, **Nitsche A** (2016). Development of a genus-specific antigen capture ELISA for orthopoxviruses – target selection and optimized screening. *PLoS One* 11: e0150110.

Hobi S, Mueller RS, Hill M, **Nitsche A** et al. (2015). Neurogenic inflammation and colliquative lymphadenitis with persistent orthopox virus DNA detection in a human case of cowpox virus infection transmitted by a domestic cat. *Br J Dermatol* 173(2): 535-539.

Medaglia ML, Moussatché N, **Nitsche A**, **Dabrowski PW**, Li Y, Damon IK, Lucas CG, Arruda LB, Damaso CR (2015). Genomic analysis, phenotype, and virulence of the historical Brazilian smallpox vaccine strain IOC: Implications for the origins and evolutionary relationships of vaccinia virus. *J Virol* 89(23): 11909-11925.

Świtaj K, Kajfasz P, **Kurth A**, **Nitsche A** (2015). Cowpox after a cat scratch - case report from Poland. *Ann Agric Environ Med* 22(3):456-458.

**Radonić A**, **Metzger S**, **Dabrowski PW**, Couacy-Hymann E, **Schuenadel L**, **Kurth A**, Mätz-Rensing K, Boesch C, Leendertz FH, **Nitsche A** (2014). Fatal monkeypox in wild-living sooty mangabey, Côte d'Ivoire, 2012. *Emerg Infect Dis* 20(6):1009-1011.

Dabrowski PW, Radonić A, **Kurth A**, **Nitsche A** (2013). Genome-wide comparison of cowpox viruses reveals a new clade related to Variola virus. *PLoS One* 8(12):e79953.

Hobi S, Mueller RS, Hill M, **Nitsche A**, Löscher T, Guggemos W, Ständer S, Rjosk-Dendorfer D, Wollenberg A (2013). Neurogenic inflammation and colliquative lymphadenitis with persistent orthopox virus DNA detection in a human case of cowpox virus infection transmitted by a domestic cat. *Br J Dermatol* 173(2):535-539.

Graef S, **Kurth A**, Auw-Haedrich C, Plange N, Kern WV, **Nitsche A**, Reinhard T (2013). Clinicopathological findings in persistent corneal cowpox infection. *JAMA Ophthalmol* 131(8):1089-1091.

Schwarzer H, **Kurth A**, Hermel M, Plange N (2013). Severe ulcerative keratitis in ocular cowpox infection. *Graefes Arch Clin Exp Ophthalmol* 251: 1451-1452.

**Kurth A**, **Nitsche A** (2012). Cowpox in Zoo Animals (chapter 5). In: R. Eric Miller, Murray E. Fowler (eds.), *Fowler's Zoo and Wild Animal Medicine: Current Therapy*, Vol. 7. St. Louis, MO: Elsevier/Saunders.

**Kurth A, Nitsche A** (2011). Detection of Human-Pathogenic Poxviruses (chapter 15). In: John R. Stephenson, Alan Warnes (eds.), *Diagnostic Virology Protocols. Methods in Molecular Biology* 665, 2. Aufl. Totowa, NJ: Humana Press.

Schupp CJ, **Nitsche A**, Bock-Hensley O, Böhm S, Flechtenmacher C, **Kurth A**, Saenger K, Hoferer M, Küsters U, Günther P, Engelmann G, Schnitzler P (2011). A 14-year-old girl with a vesicle on her finger and lymphadenitis. *J Clin Virol* 50: 1-3.

Haase O, Moser A, Rose C, **Kurth A**, Zillikens D, Schmidt E (2011). Generalised cowpox infection in a patient with Darier disease. *Br J Dermatol* 164: 1116-1118.

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**Langhammer S, Koban R, Yue C, Ellerbrok H** (2011). Inhibition of poxvirus spreading by the anti-tumor drug Gefitinib (Iressa™). *Antiviral Res* 89(1): 64-70.

**Kramski M**, Mätz-Rensing K, Stahl-Hennig C, Kaup FJ, **Nitsche A, Pauli G, Ellerbrok H** (2010). A Novel Highly Reproducible and Lethal Nonhuman Primate Model for Orthopox Virus Infection. *PLoS One* 5(4): e10412.

**Schroeder K, Nitsche A** (2010). Multicolour, multiplex real-time PCR assay for the detection of human-pathogenic poxviruses. Short Communication. *Mol Cell Probes* 24(2): 110-113.

**Radonic A, Thulke S, Achenbach J, Kurth A, Vreemann A, König T, Walter C, Possinger K, Nitsche A** (2010). Anionic polysaccharides from phototrophic microorganisms exhibit antiviral activities to vaccinia virus. *J Antivir Antiretrovir* 2: 051-055.

**Witkowski PT, Schuenadel L, Wiethaus J, Bourquain DR, Kurth A, Nitsche A** (2010). Cellular impedance measurement as a new tool for poxvirus titration, antibody neutralization testing and evaluation of antiviral substances. *Biochem Biophys Res Commun* 401(1): 37-41.

**Kurth A, Straube M, Kuczka A, Dunsche AJ, Meyer H, Nitsche A** (2009). Cowpox Virus Outbreak in Banded Mongooses (*Mungos mungo*) and Jaguarundis (*Herpailurus yagouaroundi*) with a Time-Delayed Infection to Humans. *PLoS One* 4(9): e6883.

Strenger V, Müller M, Richter S, Revilla-Fernandez S, **Nitsche A, Klee SR, Ellerbrok H, Zenz W** (2009). A 17-year-old girl with a black eschar. Cowpox virus infection. *Clin Infect Dis* 48: 91-92, 133-134.

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Becker C, **Kurth A**, Hessler F, Kramp H, Gokel M, Hoffmann R, Kuczka A, **Nitsche A** (2009). Cowpox Virus Infection in Pet Rat Owners: Not Always Immediately Recognized. Dtsch Arztebl Int 106: 329-334.

**Kurth A**, Wibbelt G, Gerber H-P, Petschaelis A, **Pauli G**, **Nitsche A** (2008). Rat-to-elephant-to-human transmission of cowpox virus [letter]. Emerg Infect Dis 14: 670-671.

Bonnekoh B, Falk K, Reckling KF, Kenklies S, **Nitsche A**, Ghebremedhin B, Pokrywka A, Franke I, Thriene B, König W, **Pauli G**, Gollnick H (2008). Cowpox infection transmitted from a domestic cat. J Dtsch Dermatol Ges 6: 210-213.

Zhang XY, **Kurth A**, Pauly D, **Pauli G**, Schade R, **Ellerbrok H** (2008). Application of hightitred IgY antibodies in orthopox virus diagnostics. J Chin Pharmaceut Sci 17(3): 183-191.

**Nitsche A**, Gelderblom HR, Eisendle K, Romani N, **Pauli G** (2007). Pitfalls in diagnosing human poxvirus infections. J Clin Virol 38: 165-168.

**Nitsche A**, **Kurth A**, **Pauli G** (2007). Viremia in human Cowpox virus infection. J Clin Virol 40(2): 160-162.

**Nitsche A**, **Pauli G** (2007). Sporadic human cases of cowpox in Germany. Euro Surv 12(4): E070419.3.

**Nitsche A**, Büttner M, Wilhelm S, **Pauli G**, Meyer H (2006). Real-time PCR detection of Parapoxvirus DNA. Clin Chem 52: 316-319.

**Nitsche A**, Steger B, **Ellerbrok H**, **Pauli G** (2005). Detection of vaccinia virus DNA on the LightCycler by fluorescence melting curve analysis. J Virol Methods 126: 187-195.

**Nitsche A**, **Ellerbrok H**, **Pauli G** (2004). Detection of orthopoxvirus DNA by real-time PCR and identification of variola virus DNA by melting analysis. J Clin Microbiol 42: 1207-1213.