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Literatur

1. AABB (2015) Listing of countries in which pathogen reduction technology systems and products are in use. <https://www.aabb.org/tm/eid/Documents/prt-systems-in-use-country-listing.pdf>. Zugriffen: 18. Mai 2018
2. Salunkhe V, van der Meer PF, de Korte D, Seghatchian J, Gutiérrez L (2015) Development of blood transfusion product pathogen reduction treatments: A review of methods, current applications and demands. *Transfus Apher Sci* 52:19–34
3. Marks DC, Faddy HM, Johnson L (2014) Pathogen reduction technologies. *ISBT Sci Ser* 9:44–50
4. Custer B (2013) Update on pathogen reduction technology. *ISBT Sci Ser* 8:80–84
5. Schlenke P (2014) Pathogen inactivation technologies for cellular blood components: an update. *Transfus Med Hemother* 41:309–325
6. REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 5 April 2017 on medical devices, amending Directive 2001/83/EC, Regulation (EC) No 178/2002 and Regulation (EC) No 1223/2009 and repealing Council Directives 90/385/EEC and 93/42/EEC

- <http://data.europa.eu/eli/reg/2017/745/oj>
Zugegriffen: 18. 05. 2018
7. Funk MB, Lohmann A, Spranger R (2015) Hämovigilanz-Bericht des Paul-Ehrlich-Instituts 2013/14. <http://www.pei.de/haemovigilanzbericht>. Zugegriffen: 18. Mai 2018
 8. Funk MB, Heiden M (2017) Hämovigilanz-Bericht des Paul-Ehrlich-Instituts 2015: Auswertung der Meldungen von schwerwiegenden Reaktionen und Zwischenfällen nach § 63i AMG. <http://www.pei.de/haemovigilanzbericht>. Zugegriffen: 18. Mai 2018
 9. Walther-Wenke G, Wirsing von König CH, Däubener W et al (2011) Monitoring bacterial contamination of blood components in Germany: effect of contamination reduction measures. *Vox Sang* 100:359–366
 10. Irsch J, Lin L (2011) Pathogen inactivation of platelet and plasma blood components for transfusion using the INTERCEPT™ Blood System. *Transfus Med Hemother* 38:19–31
 11. CERUS Pressemitteilung vom 16.12.2014. <http://www.cerus.com/Investors/Press-Releases/Press-Release-Details/2014/FDA-Approves-INTERCEPT-Blood-System-for-Plasma/default.aspx>. Zugegriffen: 18. Mai 2018
 12. CERUS Pressemitteilung vom 18.12.2014. <http://www.cerus.com/Investors/Press-Releases/Press-Release-Details/2014/FDA-Approves-INTERCEPT-Blood-System-for-Platelets/default.aspx>. Zugegriffen: 18. Mai 2018
 13. Liu W, Cimino GD, Corash L, Lin L (2011) The extent of amotosalen photodegradation during photochemical treatment of platelet components correlates with the level of pathogen inactivation. *Transfusion* 51:52–61
 14. Hanson CV (1992) Photochemical inactivation of viruses with psoralens: an overview. *Blood Cells* 18:7–25
 15. Wollowitz S (2001) Fundamentals of the psoralen-based Helinx technology for inactivation of infectious pathogens and leukocytes in platelets and plasma. *Semin Hematol* 38:4–11
 16. Kanne D, Straub K, Hearst JE, Rapoport H (1982) Isolation and characterization of pyrimidine-psoralen-pyrimidine photodiadducts from DNA. *J Am Chem Soc* 104:6754–6764
 17. Prudent M, Sonogo G, Abonnenc M, Tissot JD, Lion N (2014) LC-MS/MS Analysis and Comparison of Oxidative Damages on Peptides Induced by Pathogen Reduction Technologies for Platelets. *J Am Soc Mass Spectrom* 25:651–661
 18. Sonogo G, Abonnenc M, Tissot JD, Prudent M, Lion N (2017) Redox Proteomics and Platelet Activation: Understanding the Redox Proteome to Improve Platelet Quality for Transfusion. *Int J Mol Sci* 18:387–408
 19. Lin L, Londe H, Janda JM, Hanson CV, Corash L (1994) Photochemical inactivation of pathogenic bacteria in human platelet concentrates. *Blood* 83:2698–2706
 20. Lin L, Cook ND, Wiesehahn GP et al (1997) Photochemical inactivation of viruses and bacteria in platelet concentrates by use of a novel psoralen and long-wavelength ultraviolet light. *Transfusion* 37:423–435
 21. Lin L (2001) Inactivation of cytomegalovirus in platelet concentrates using Helinx technology. *Semin Hematol* 38:27–33
 22. Lin L, Dikeman R, Molini B et al (2004) Photochemical treatment of platelet concentrates with amotosalen and long-wavelength ultraviolet light inactivates a broad spectrum of pathogenic bacteria. *Transfusion* 44:1496–1504
 23. Lin L, Hanson CV, Alter HJ et al (2005) Inactivation of viruses in platelet concentrates by photochemical treatment with amotosalen and long-wavelength ultraviolet light. *Transfusion* 45:580–590
 24. Eastman RT, Barrett LK, Dupuis K, Buckner FS, Van Voorhis WC (2005) Leishmania inactivation in human pheresis platelets by a psoralen (amotosalen HCl) and long-wavelength ultraviolet irradiation. *Transfusion* 45:1459–1463
 25. Jauvin V, Alfonso RD, Guillemin B, Dupuis K, Fleury HJ (2005) In vitro photochemical inactivation of cell-associated human T-cell leukemia virus Type I and II in human platelet concentrates and plasma by use of amotosalen. *Transfusion* 45:1151–1159
 26. Pinna D, Sampson-Johannes A, Clementi M et al (2005) Amotosalen photochemical inactivation of severe acute respiratory syndrome coronavirus in human platelet concentrates. *Transfus Med* 15:269–276
 27. Gallian P, Vignoli C, Dombey AM et al (2006) Inactivation of a European strain of West Nile virus in single-donor platelet concentrate using the INTERCEPT blood system. *Vox Sang* 91:345–347
 28. Roback JD, Conlan M, Drew WL, Ljungman P, Nichols WG, Preiksaitis JK (2006) The role of photochemical treatment with amotosalen and UV-A light in the prevention of transfusion-transmitted cytomegalovirus infections. *Transfus Med Rev* 20:45–56
 29. Brecher ME, Hay S, Corash L, Hsu J, Lin L (2007) Evaluation of bacterial inactivation in prestorage pooled, leukoreduced, whole blood-derived platelet concentrates suspended in plasma prepared with photochemical treatment. *Transfusion* 47:1896–1901
 30. Nussbaumer W, Allerstorfer D, Grabmer C et al (2007) Prevention of transfusion of platelet components contaminated with low levels of bacteria: a comparison of bacteria culture and pathogen inactivation methods. *Transfusion* 47:1125–1133
 31. Grellier P, Benach J, Labaied M et al (2008) Photochemical inactivation with amotosalen and long-wavelength ultraviolet light of Plasmodium and Babesia in platelet and plasma components. *Transfusion* 48:1676–1684
 32. Kwon SY, Kim IS, Bae JE et al (2014) Pathogen inactivation efficacy of MIRASOL PRT System and INTERCEPT Blood System for non-leukoreduced platelet-rich plasma-derived platelets suspended in plasma. *Vox Sang* 107:254–260
 33. Musso D, Richard V, Broult J, Cao-Lormeau VM (2014) Inactivation of dengue virus in plasma with amotosalen and ultraviolet A illumination. *Transfusion* 54:2924–2930
 34. Aubry M, Richard V, Green J, Broult J, Musso D (2016) Inactivation of Zika virus in plasma with amotosalen und ultraviolet A illumination. *Transfusion* 56:33–40
 35. CERUS (2016) INTERCEPT™ Plasma Technical Data Sheet. http://interceptbloodsystem.com/sites/default/files/resources/technical_data_sheet_-_plasma_-_english.pdf. Zugegriffen: 18. Mai 2018
 36. CERUS (2017) INTERCEPT™ Platelets Technical Data Sheet. http://interceptbloodsystem.com/sites/default/files/resources/prd-tds_00121_v9.0_secure_0.pdf. Zugegriffen: 18. Mai 2018
 37. Hauser L, Roque-Alfonso AM, Beylouné A et al (2014) Hepatitis E transmission by transfusion of INTERCEPT blood system-treated plasma. *Blood* 123:796–797
 38. Juhl D, Baylis SA, Blümel J, Görg S, Hennig H (2014) Seroprevalence and incidence of hepatitis E virus infection in German blood donors. *Transfusion* 54:49–56
 39. Panning M, Umhau M, Emmerich F, Wedemeyer H, Pischke S (2015) Diskussion zu dem Beitrag Hepatitis E in Deutschland-eine unterschätzte Infektionskrankheit. *Dtsch Arztebl* 112:220
 40. Baylis SA, Corman VM, Ong E, Linnen JM, Nübling CM, Blümel J (2016) Hepatitis E viral loads in plasma pools for fractionation. *Transfusion* 56:2532–2537
 41. Domanovic D, Tedder R, Blümel J et al (2017) Hepatitis E and blood donation safety in selected European countries: a shift to screening? *EuroSurveillance* 22:30514
 42. Schmidt M, Hourfar MK, Sireis W et al (2015) Evaluation of the effectiveness of a pathogen inactivation technology against clinically relevant transfusion-transmitted bacterial strains. *Transfusion* 55:2104–2112
 43. Störmer M, Arroyo A, Brachert J et al (2012) Establishment of the first international repository for transfusion-relevant bacteria reference strains: ISBT working party transfusion-transmitted infectious diseases (WP-TTID), subgroup on bacteria. *Vox Sang* 102:22–31
 44. Benjamin RJ, Wagner SJ (2015) Bacterial pathogen reduction requires validation under conditions of intended use. *Transfusion* 55:2060–2063
 45. Grass JA, Hei DJ, Metchette K et al (1998) Inactivation of leukocytes in platelet concentrates by photochemical treatment with psoralen plus UVA. *Blood* 91:2180–2188
 46. Grass JA, Wafa T, Reames A et al (1999) Prevention of transfusion-associated graft-versus-host disease by photochemical treatment. *Blood* 93:3140–3147
 47. Tice RR, Gatehouse D, Kirkland D, Speit G (2007) The pathogen reduction treatment of platelets with S-59 HCl (Amotosalen) plus ultraviolet A light: Genotoxicity profile and hazard assessment. *Mut Res* 630:50–68
 48. Lin L, Conlan MG, Tessman J, Cimino G, Porter S (2005) Amotosalen interactions with platelet and plasma components: absence of neoantigen formation after photochemical treatment. *Transfusion* 45:1610–1620
 49. Van Rhenen DJ, Vermeij J, Mayaudon V, Hind C, Lin L, Corash L (2000) Functional characteristics of S-59 photochemically treated platelet concentrates derived from buffy coats. *Vox Sang* 79:206–214
 50. Jansen GA, van Vlieth HH, Vermeij H et al (2004) Functional characteristics of photochemically treated platelets. *Transfusion* 44:313–319
 51. Van Aelst B, Feys HB, Devloo R, Vanhoorelbeke K, Vandekerckhove P, Compennolle V (2015) Riboflavin and amotosalen photochemical treatments of platelet concentrates reduce thrombus formation kinetics in vitro. *Vox Sang* 108:328–339
 52. Apelseth TO, Bruserud O, Wentzel-Larsen T, Bakken AM, Bjorsvik S, Hervig T (2007) In vitro evaluation of metabolic changes and residual platelet responsiveness in photochemical treated and gamma-irradiated single-donor platelet concentrates during long-term storage. *Transfusion* 47:653–665

53. Johnson L, Loh YS, Kwok M, Marks DS (2013) In vitro assessment of buffy-coat derived platelet components suspended in SSP+ treated with the INTERCEPT Blood system. *Transfus Med* 23:121–129
54. Janetzko K, Lin L, Eichler H, Mayaudon V, Flament J, Klüter H (2004) Implementation of the INTERCEPT Blood System for Platelets into routine blood bank manufacturing procedures: evaluation of apheresis platelets. *Vox Sang* 86:239–245
55. Picker SM, Speer R, Gathof BS (2004) Functional characteristics of buffy-coat PLTs photochemically treated with amotosalen-HCl for pathogen inactivation. *Transfusion* 44:320–329
56. Sandgren P, Diedrich B (2015) Pathogen inactivation of double-dose buffy-coat platelet concentrates photochemically treated with amotosalen and UVA light: preservation of in vitro function. *Vox Sang* 108:340–349
57. Lozano M, Galan A, Mazzara R, Corash L, Escolar G (2007) Leukoreduced buffy coat-derived platelet concentrates photochemically treated with amotosalen HCl and ultraviolet A light stored up to 7 days: assessment of hemostatic function under flow conditions. *Transfusion* 47:666–671
58. Picker SM, Schneider V, Gathof BS (2009) Platelet function assessed by shear-induced deposition of split triple-dose apheresis concentrates treated with pathogen reduction technologies. *Transfusion* 49:1224–1232
59. Hechler B, Ohlmann P, Chafey P et al (2013) Preserved functional and biochemical characteristics of platelet components prepared with amotosalen and ultraviolet A for pathogen inactivation. *Transfusion* 53:1187–1200
60. Tynngard N, Johansson BM, Lindahl TL, Berlin G, Hansson M (2008) Effects of INTERCEPT pathogen inactivation on platelet function as analysed by free oscillation rheometry. *Transfus Apher Sci* 38:85–88
61. Osman A, Hitzler WE, Meyer CU et al (2015) Effects of pathogen reduction systems on platelet microRNAs, mRNAs, activation, and function. *Platelets* 26:154–163
62. Bakkour S, Chafets DM, Wen L et al (2016) Assessment of nucleic acid modification induced by amotosalen and ultraviolet A light treatment of platelets and plasma using real-time polymerase chain reaction amplification of variable length fragments of mitochondrial DNA. *Transfusion* 56:410–420
63. Laffont B, Corduan A, Plé H et al (2013) Activated platelets can deliver mRNA regulatory Ago2-microRNA complexes to endothelial cells via microparticles. *Blood* 122:253–261
64. Semple JW (2013) Platelets deliver small packages of genetic function. *Blood* 122:155–156
65. Osman A, Hitzler WE, Ameur A, Provost P (2015) Differential Expression Analysis by RNA-Seq Reveals Perturbations in the Platelet mRNA Transcriptome Triggered by Pathogen Reduction Systems. *Plos One* 10:e133070
66. Prudent M, Crettaz D, Delobel J, Tissot JD, Lion N (2012) Proteomic analysis of Intercept-treated platelets. *J Proteomics* 76:316–328
67. Prudent M, D'Alessandro A, Cazenave JP et al (2014) Proteome changes in platelets after pathogen inactivation—an interlaboratory consensus. *Transfus Med Rev* 28:72–83
68. Terumo BCT MIRASOL® Pathogen Reduction Technology System, 05.12.2014. <https://www.terumobct.com/mirasol>. Zugegriffen: 22. Mai 2018
69. Corbin F (2002) Pathogen inactivation of blood components: current status and introduction of an approach using riboflavin as a photosensitizer. *Int J Hematol* 76(Suppl 2):253–257
70. Marschner S, Goodrich R (2011) Pathogen reduction technology treatment of platelets, plasma and whole blood using riboflavin and UV light. *Transfus Med Hemother* 38:8–18
71. Owada T, Kaneko M, Matsumoto C et al (2014) Establishment of culture systems for genotypes 3 and 4 hepatitis E virus (HEV) obtained from human blood and application of HEV inactivation using a pathogen reduction technology system. *Transfusion* 54:2820–2827
72. Ruane PH, Edrich R, Gamp D, Keil SD, Leonard RL, Goodrich RP (2004) Photochemical inactivation of selected viruses and bacteria in platelet concentrates using riboflavin and light. *Transfusion* 44:877–885
73. Keil SD, Bengrine A, Bowen R et al (2015) Inactivation of viruses in platelet and plasma products using a riboflavin-and-UV-based photochemical treatment. *Transfusion* 55:1736–1744
74. Faddy HM, Fryk JJ, Watterson D et al (2016) Riboflavin and ultraviolet light: impact on dengue virus infectivity. *Vox Sang* 111:235–241
75. Faddy HM, Prow NA, Fryk JJ et al (2015) The effect of riboflavin and ultraviolet light on the infectivity of arboviruses. *Transfusion* 55:824–831
76. Vanlandingham DL, Keil SD, Horne KM, Pyles R, Goodrich RP, Higgs S (2013) Photochemical inactivation of chikungunya virus in plasma and platelets using the MIRASOL pathogen reduction technology system. *Transfusion* 53:284–290
77. Keil SD, Saakadze N, Bowen R et al (2015) Riboflavin and ultraviolet light for pathogen reduction of murine cytomegalovirus in blood products. *Transfusion* 55:858–863
78. Goodrich RP, Edrich RA, Li J, Seghatchian J (2006) The MIRASOL™ PRT system for pathogen reduction of platelets and plasma: an overview of current status and future trends. *Transfus Apher Sci* 35:5–17
79. Shinohara N, Matsumoto C, Chatani M et al (2015) Efficacy of the MIRASOL pathogen reduction technology system against severe fever with thrombocytopenia syndrome virus (SFTSV). *Vox Sang* 109:417–419
80. Rentas F, Harman R, Gomez C et al (2007) Inactivation of *Orientia tsutsugamushi* in red blood cells, plasma, and platelets with riboflavin and light, as demonstrated in an animal model. *Transfusion* 47:240–247
81. Tonnetti L, Proctor MC, Reddy HL, Goodrich RP, Leiby DA (2010) Evaluation of the MIRASOL pathogen reduction technology system against *Babesia microti* in apheresis platelets and plasma. *Transfusion* 50:1019–1027
82. Tonnetti L, Thorp AM, Reddy HL, Keil SD, Goodrich RP, Leiby DA (2013) Riboflavin and ultraviolet light reduce the infectivity of *Babesia microti* in whole blood. *Transfusion* 53:860–867
83. Tonnetti L, Thorp AM, Reddy HL et al (2015) Reduction of *Leishmania donovani* infectivity in whole blood using riboflavin and ultraviolet light. *Transfusion* 55:326–329
84. El Chaar M, Atwal S, Freimanis GL, Dinko B, Sutherland CJ, Allain JP (2013) Inactivation of *Plasmodium falciparum* in whole blood by riboflavin plus irradiation. *Transfusion* 53:3174–3183
85. Keil SD, Kiser P, Sullivan JJ et al (2013) Inactivation of *Plasmodium* spp. in plasma and platelet concentrates using riboflavin and ultraviolet light. *Transfusion* 53:2278–2286
86. Tonnetti L, Thorp AM, Reddy HL, Keil SD, Goodrich RP, Leiby DA (2012) Evaluating pathogen reduction of *Trypanosoma cruzi* with riboflavin and ultraviolet light for whole blood. *Transfusion* 52:409–416
87. Cardo LJ, Salata J, Mendez J, Reddy H, Goodrich R (2007) Pathogen inactivation of *Trypanosoma cruzi* in plasma and platelet concentrates using riboflavin and ultraviolet light. *Transfus Apher Sci* 37:131–137
88. Brecher ME, Hay SN (2005) Bacterial contamination of blood components. *Clin Microbiol Rev* 18:195–204
89. Goodrich RP, Gilmour D, Hovenga N, Keil SD (2009) A laboratory comparison of pathogen reduction technology treatment and culture of platelet products for addressing bacterial contamination concerns. *Transfusion* 49:1205–1216
90. Fast LD, DiLeone G, Li J, Goodrich R (2006) Functional inactivation of white blood cells by MIRASOL treatment. *Transfusion* 46:642–648
91. Fast LD, DiLeone G, Cardarelli G, Li J, Goodrich R (2006) MIRASOL PRT treatment of donor white blood cells prevents the development of xenogeneic graft-versus-host disease in Rag2^{-/-}γc^{-/-} double knockout mice. *Transfusion* 46:1553–1560
92. Jackman RP, Heitman JW, Marschner S, Goodrich RP, Norris PJ (2009) Understanding loss of donor white blood cell immunogenicity after pathogen reduction: mechanisms of action in ultraviolet illumination and riboflavin treatment. *Transfusion* 49:2686–2699
93. Marschner S, Fast LD, Baldwin WM III, Slichter SJ, Goodrich RP (2010) White blood cell inactivation after treatment with riboflavin and ultraviolet light. *Transfusion* 50:2489–2498
94. Fast LD, DiLeone G, Marschner S (2011) Inactivation of human white blood cells in platelet products after pathogen reduction technology treatment in comparison to gamma irradiation. *Transfusion* 51:1397–1404
95. Fast LD, Nevola M, Tavares J, Reddy HL, Goodrich RP, Marschner S (2013) Treatment of whole blood with riboflavin plus ultraviolet light, an alternative to gamma irradiation in the prevention of transfusion-associated graft-versus-host disease? *Transfusion* 53:373–381
96. Reddy HL, Doane SK, Keil SD, Marschner S, Goodrich RP (2013) Development of a riboflavin and ultraviolet light-based device to treat whole blood. *Transfusion* 53:1315–1365
97. FDA Database of Select Committee on GRAS Substances (SCOGS) Reviews: Riboflavin. <https://www.fda.gov/Food/IngredientsPackagingLabeling/GRAS/SCOGS/default.htm>. Zugegriffen: 22. Mai 2018
98. Diakonis VF, Grentzelos MA, Tzatzarakis MN et al (2012) Riboflavin's time-dependent degradation rate induced by ultraviolet A irradiation. *Eur J Ophthalmol* 22(Suppl. 7):51–56
99. Ambruso DR, Thurman G, Marschner S, Goodrich RP (2009) Lack of antibody formation to platelet neoantigens after transfusion of riboflavin and ultraviolet light-treated platelet concentrates. *Transfusion* 49:2631–2636
100. Janetzko K, Hinz K, Marschner S, Goodrich R, Klüter H (2009) Pathogen reduction technology (MIRASOL) treated single-donor platelets

- resuspended in a mixture of autologous plasma and PAS. *Vox Sang* 97:234–239
101. Li J, Goodrich L, Hansen E, Edrich R, Gampp D, Goodrich RP (2005) Platelet glycolytic flux increases stimulated by ultraviolet-induced stress is not the direct cause of platelet morphology and activation changes: possible implications for the role of glucose in platelet storage. *Transfusion* 45:1750–1758
 102. Picker SM, Tauszig ME, Gathof BS (2012) Cell quality of apheresis-derived platelets treated with riboflavin-ultraviolet light after resuspension in platelet additive solution. *Transfusion* 52:510–516
 103. Ostrowski SR, Bochsén L, Salado-Jimena JA et al (2010) In vitro cell quality of buffy coat platelets in additive solution treated with pathogen reduction technology. *Transfusion* 50:2210–2219
 104. Johnson L, Winter KM, Reid S et al (2011) The effect of pathogen reduction technology (MIRASOL) on platelet quality when treated in additive solution with low plasma carryover. *Vox Sang* 101:208–214
 105. Cookson P, Thomas S, Marschner S, Goodrich R, Cardigan R (2012) In vitro quality of single-donor platelets treated with riboflavin and ultraviolet light and stored in platelet storage medium for up to 8 days. *Transfusion* 52:983–994
 106. Perez-Pujol S, Tonda R, Lozano M et al (2005) Effects of a new pathogen-reduction technology (MIRASOL PRT) on functional aspects of platelet concentrates. *Transfusion* 45:911–919
 107. Picker SM, Steisel A, Gathof BS (2008) Effects of MIRASOL PRT treatment on storage lesion development in plasma-stored apheresis-derived platelets compared to untreated and irradiated units. *Transfusion* 48:1685–1692
 108. Picker SM, Steisel A, Gathof BS (2009) Evaluation of White Blood Cell- and Platelet-Derived Cytokine Accumulation in MIRASOL-PRT-Treated Platelets. *Transfus Med Hemother* 36:114–120
 109. Picker SM, Oustianskaia L, Schneider V, Gathof BS (2010) Annexin V Release and Transmembrane Mitochondrial Potential during Storage of Apheresis-Derived Platelets Treated for Pathogen Reduction. *Transfus Med Hemother* 37:7–12
 110. Castrillo A, Cardoso M, Rouse L (2013) Treatment of Buffy Coat Platelets in Platelet Additive Solution with the MIRASOL® Pathogen Reduction Technology System. *Transfus Med Hemother* 40:44–48
 111. Zeddis S, De Cuyper IM, van der Meer PF et al (2014) Pathogen reduction treatment using riboflavin and ultraviolet light impairs platelet reactivity toward specific agonists in vitro. *Transfusion* 54:2292–2300
 112. Galan AM, Lozano M, Molina P et al (2011) Impact of pathogen reduction technology and storage in platelet additive solutions on platelet function. *Transfusion* 51:808–815
 113. Terada C, Mori J, Okazaki H, Satake M, Tadokoro K (2014) Effects of riboflavin and ultraviolet light treatment on platelet thrombus formation on collagen via integrin α IIb β 3 activation. *Transfusion* 54:1808–1816
 114. Klein-Bosgoed C, Schubert P, Devine DV (2016) Riboflavin and ultraviolet illumination affects selected platelet mRNA transcript amounts differently. *Transfusion* 56:2286–2295
 115. Feys HB, van Aelst B, Devreese K et al (2014) Oxygen removal during pathogen inactivation with riboflavin and UV light preserves protein function in plasma for transfusion. *Vox Sang* 106:307–315
 116. Johnson L, Marks D (2015) Treatment of Platelet Concentrates with the MIRASOL Pathogen Inactivation System Modulates Platelet Oxidative Stress and NF- κ B Activation. *Transfus Med Hemother* 42:167–173
 117. Schubert P, Culibrk B, Coupland D, Scammell K, Gyongyossy-Issa M, Devine DV (2012) Riboflavin and ultraviolet light treatment potentiates vasodilator-stimulated phosphoprotein Ser-239 phosphorylation in platelet concentrates during storage. *Transfusion* 52:397–408
 118. AuBuchon JP, Herschel L, Roger J et al (2005) Efficacy of apheresis platelets treated with riboflavin and ultraviolet light for pathogen reduction. *Transfusion* 45:1335–1341
 119. Goodrich RP, Li J, Pieters H, Crookes R, Roodt J, Heyns Adu P (2006) Correlation of in vitro platelet quality measurements with in vivo platelet viability in human subjects. *Vox Sang* 90:279–285
 120. Macopharma (2014) Macopharma Theraflex UV-Platelets. <http://blood-safety.macopharma.com/en/category/documents-literatures/theraflex-uv-platelets>. Zugegriffen: 22. Mai 2018
 121. Seltam A, Müller TH (2011) UVC Irradiation for Pathogen Reduction of Platelet Concentrates and Plasma. *Transfus Med Hemother* 38:43–54
 122. Seghatchian J, Tolksdorf F (2012) Characteristics of the THERAFLEX UV-Platelets pathogen inactivation system - an update. *Transfus Apher Sci* 46:221–229
 123. Mohr H, Steil L, Gravemann U et al (2009) A novel approach to pathogen reduction in platelet concentrates using short-wave ultraviolet light. *Transfusion* 49:2612–2624
 124. Thiele T, Pohler P, Kohlmann T et al (2015) Tolerance of platelet concentrates treated with UVC-light only for pathogen reduction – a phase I clinical trial. *Vox Sang* 109:44–51
 125. Kim S, Handke W, Gravemann U et al (2018) Mitochondrial DNA multiplex real-time polymerase chain reaction inhibition assay for quality control of pathogen inactivation by ultraviolet C light in platelet concentrates. *Transfusion* 58:758–765
 126. Steinmann E, Gravemann U, Friesland M et al (2013) Two pathogen reduction technologies – methylene blue plus light and shortwave ultraviolet light – effectively inactivate hepatitis C virus in blood products. *Transfusion* 53:1010–1018
 127. van der Meer PF, Gravemann U, de Korte D et al (2016) Effect of increased agitation speed on pathogen inactivation efficacy and in vitro quality in UVC-treated platelet concentrates. *Vox Sang* 111:127–134
 128. Faddy HM, Fryk JJ, Prow NA et al (2016) Inactivation of dengue, chikungunya, and Ross River viruses in platelet concentrates after treatment with ultraviolet C light. *Transfusion* 56:1548–1555
 129. Seltam A (2011) Pathogenreduktion mittels UVC-Licht: Entwicklung des THERAFLEX UV-Plättchen-Systems. DGTI, Freiburg. <https://www.dgti.de/docs/doclink/10802/6.%20Vorortrag%20Seltam,%20Theraflex%20UV%20Platelets.pdf>. Zugegriffen: 23. Mai 2018
 130. Gravemann U, Schmidt JP, Tolksdorf F et al (2013) Influenza A virus H3N2 is efficiently inactivated by the THERAFLEX UV-PLATELETS system. *Vox Sang* 105(Suppl. 1):134
 131. Fryk JJ, Marks DC, Hobson-Peters J et al (2017) Reduction of Zika virus infectivity in platelet concentrates after treatment with ultraviolet C light and in plasma after treatment with methylene blue and visible light. *Transfusion* 57:2677–2682
 132. Gravemann U, Lambrecht B, Schmidt J-P, Seltam A (2014) Hepatitis A virus is efficiently inactivated by the THERAFLEX UV-Platelets system. *Vox Sang* 107(Suppl. 1):128
 133. Spindler-Raffel E, Benjamin RJ, McDonald CP et al (2017) ISBT Working Party Transfusion-Transmitted Infectious Diseases (WP-TTID), Subgroup on Bacteria (2017) Enlargement of the WHO international repository for platelet transfusion-relevant bacteria reference strains. *Vox Sang* 112:713–722
 134. Gravemann U, Tolksdorf F, Handke W, Müller TH, Seltam A (2017) The THERAFLEX UV-Platelets technology efficiently inactivates transfusion-relevant bacteria species in contaminated platelet concentrates. *Transfusion* 57(Suppl. 1):203
 135. Castro E, Girones N, Guerrero N, Barea L, Fresno M (2008) The effectiveness of UVC pathogen inactivation system on reducing the Trypanosoma cruzi and Leishmania infantum burden in platelets. *Vox Sang* 95(Suppl. 1):290
 136. Castro E, González LM, Bautista JM, Montero E, Rubio JM (2013) The efficacy of UVC pathogen inactivation on the reduction of Plasmodium falciparum in buffy coat derived platelets. *Vox Sang* 105(Suppl. 1):133
 137. Castro E, González LM, Rubio JM, Ramiro R, Gironés N, Montero E (2014) The efficacy of the ultraviolet C pathogen inactivation system in the reduction of Babesia divergens in pooled buffy coat platelets. *Transfusion* 54:2207–2216
 138. Pohler P, Müller M, Winkler C et al (2015) Pathogen reduction by ultraviolet C light effectively inactivates human white blood cells in platelet products. *Transfusion* 55:337–347
 139. Pohler P, Lehmann J, Veneruso V et al (2012) Evaluation of the tolerability and immunogenicity of ultraviolet C-irradiated autologous platelets in a dog model. *Transfusion* 52:2414–2426
 140. Bashir S, Cookson P, Wiltshire M et al (2013) Pathogen inactivation of platelets using ultraviolet C light: effect on in vitro function and recovery and survival of platelets. *Transfusion* 53:990–1000
 141. Tynngard N, Trinks M, Berlin G (2015) In vitro function of platelets treated with ultraviolet C light for pathogen inactivation: a comparative study with non irradiated and gamma-irradiated platelets. *Transfusion* 55:1169–1177
 142. Van Aelst B, Devloo R, Vandekerckhove P, Compennolle V, Feys HB (2015) Ultraviolet C light pathogen inactivation treatment of platelet concentrates preserves integrin activation but affects thrombus formation kinetics on collagen in vitro. *Transfusion* 55:2404–2414
 143. Martini R, Kempfer CB, Rodrigues Mde A et al (2010) Bacterial contamination on platelet concentrates: identification, antimicrobial susceptibility profile and transfusion-related sepsis. *Rev Soc Bras Med Trop* 43:682–685
 144. Kou Y, Pagotto F, Hannach B, Ramirez-Arcos S (2015) Fatal false-negative transfusion infection involving a buffy coat platelet pool contaminated with biofilm-positive Staphylococcus epidermidis: a case report. *Transfusion* 55:2384–2389
 145. Snyder E, Raife T, Lin L et al (2004) Recovery and life span of ¹¹¹Indium-radiolabeled platelets treated with pathogen inactivation with amotosalen HCl (S-59) and ultraviolet A light. *Transfusion* 44:1732–1740

146. Murphy S (2004) Radiolabeling of PLTs to assess viability: a proposal for a standard. *Transfusion* 44:131–133
147. FDA Blood Products Advisory Committee (2004) July 22–23th, Summary Topic II: New Standards for Platelet Evaluation. <https://wayback.archive-it.org/7993/20170404063056/https://www.fda.gov/ohrms/dockets/ac/04/transcripts/2004-4057t1.htm>. Zugegriffen: 23. Mai 2018
148. Swiss Medic Haemovigilance Annual report 2015. www.swissmedic.ch/haemovigilance-report. Zugegriffen: 23. Mai 2018
149. Butler C, Doree C, Estcourt LJ et al (2013) Pathogen-reduced platelets for the prevention of bleeding. *Cochrane Database Syst Rev* 3:CD9072
150. Estcourt LJ, Malouf R, Hopewell S et al (2017) Pathogen-reduced platelets for the prevention of bleeding. *Cochrane Database Syst Rev* 7:CD9072
151. McCullough J, Vesole DH, Benjamin RJ et al (2004) Therapeutic efficacy and safety of platelets treated with a photochemical process for pathogen inactivation: the SPRINT Trial. *Blood* 104:1534–1541
152. Kerkhoffs JLH, van Putten WLJ, Novotny VMJ et al (2010) Clinical effectiveness of leucoreduced, pooled donor platelet concentrates, stored in plasma or additive solution with and without pathogen reduction. *Br J Haematol* 150:209–217
153. Cid J, Escolar G, Lozano M (2012) Therapeutic efficacy of platelet components treated with amotosalen and ultraviolet A pathogen inactivation method: results of a meta-analysis of randomized controlled trials. *Vox Sang* 103:322–330
154. Garban F, Guyard A, Labussière H et al (2018) Comparison of the Hemostatic Efficacy of Pathogen-Reduced Platelets vs Untreated Platelets in Patients With Thrombocytopenia and Malignant Hematologic Diseases: A Randomized Clinical Trial. *JAMA Oncol* 4:468–475
155. The MIRASOL Clinical Evaluation Study Group (2010) A randomized controlled clinical trial evaluating the performance and safety of platelets treated with MIRASOL pathogen reduction technology. *Transfusion* 50:2362–2375
156. Snyder E, McCullough J, Slichter SJ et al (2005) Clinical safety of platelets photochemically treated with amotosalen HCl and ultraviolet A light for pathogen inactivation: the SPRINT trial. *Transfusion* 45:1864–1875
157. Kleinman S, Reed W, Stassinopoulos A (2013) A patient-oriented risk-benefit analysis of pathogen-inactivated blood components: application to apheresis platelets in the United States. *Transfusion* 53:1603–1618
158. Amato M, Schennach H, Astl M et al (2017) Impact of platelet pathogen inactivation on blood component utilization and patient safety in a large Austrian Regional Medical Centre. *Vox Sang* 112:47–55
159. Knutson F, Osselaer J, Pierelli L et al (2015) A prospective, active haemovigilance study with combined cohort analysis of 19,175 transfusions of platelet components prepared with amotosalen-UVA photochemical treatment. *Vox Sang* 109:343–352
160. Jutzi M, Rüesch M, Taleghani BM (2013) Einführung der Pathogeninaktivierung für Thrombozytenkonzentrate in der Schweiz. *Schweiz Med Forum* 13:222–226