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Referenzen:

Referenzen:

- 1 AABB: (2015) Listing of countries in which pathogen reduction technology systems and products are in use. Updated March <https://www.aabb.org/tm/eid/Documents/prt-systems-in-use-country-listing.pdf> Zugegriffen: 18.05.2018
- 2 Salunkhe V, van der Meer PF, de Korte D, Seghatchian J, Gutiérrez L (2015) Development of blood transfusion product pathogen reduction treatments: A review of methods, current applications and demands. *Transfus Apher Sci* 52:19-34
- 3 Marks DC, Faddy HM, Johnson L (2014) Pathogen reduction technologies. *ISBT Science Series* 9:44-50
- 4 Custer B (2013) Update on pathogen reduction technology. *ISBT Science Series* 8: 80-84
- 5 Schlenke P (2014) Pathogen inactivation technologies for cellular blood components: an update. *Transfus Med Hemother* 41:309-325
- 6 REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 5 April 2017 on medical devices, amending Directive 2001/83/EC, Regulation (EC) No 178/2002 and Regulation (EC) No 1223/2009 and repealing Council Directives 90/385/EEC and 93/42/EEC <http://data.europa.eu/eli/reg/2017/745/oj> Zugegriffen: 18. 05. 2018
- 7 Funk MB, Lohmann A, Spranger R. Hämovigilanz-Bericht des Paul-Ehrlich-Instituts 2013/14. <http://www.pei.de/haemovigilanzbericht> Zugegriffen 18.05.2018
- 8 Funk MB, Heiden M Hämovigilanz-Bericht des Paul-Ehrlich-Instituts 2015: Auswertung der Meldungen von schwerwiegenden Reaktionen und Zwischenfällen nach § 63i AMG. <http://www.pei.de/haemovigilanzbericht> Zugegriffen 18.05.2018
- 9 Walther-Wenke G, Wirsing von König CH, Däubener W et al. (2011) Monitoring bacterial contamination of blood components in Germany: effect of contamination reduction measures. *Vox Sang* 100:359-366
- 10 Irsch J, Lin L (2011) Pathogen inactivation of platelet and plasma blood components for transfusion using the INTERCEPT™ Blood System. *Trans Med Hemother* 38:19-31
- 11 CERUS, Pressemitteilung vom 16.12.2014 <http://www.cerus.com/Investors/Press-Releases/Press-Release-Details/2014/FDA-Approves-INTERCEPT™-Blood-System-for-Plasma/default.aspx> Zugegriffen 18.05.2018
- 12 CERUS, Pressemitteilung vom 18.12.2014 <http://www.cerus.com/Investors/Press-Releases/Press-Release-Details/2014/FDA-Approves-INTERCEPT™-Blood-System-for-Platelets/default.aspx> Zugegriffen 18.05.2018
- 13 Liu W, Cimino GD, Corash L, Lin L (2011) The extent of amotosalen photodegradation during photochemical treatment of platelet components correlates with the level of pathogen inactivation. *Transfusion* 51:52-61
- 14 Hanson CV (1992) Photochemical inactivation of viruses with psoralens: an overview. *Blood Cells* 18:7-25.
- 15 Wollowitz S (2001) Fundamentals of the psoralen-based Helinx technology for inactivation of infectious pathogens and leukocytes in platelets and plasma. *Semin Hematol* 38:4-11
- 16 Kanne D, Straub K Hearst JE at al. (1982) Isolation and characterization of pyrimidine-psoralen, pyrimidine photodiadducts from DNA. *J Am Chem Soc* 104:6754-6764
- 17 Prudent M, Sonogo G, Abonnenc M, Tissot J-D, Lion N (2014) LC-MS/MS Analysis and Comparison of Oxidative Damages on Peptides Induced by Pathogen Reduction Technologies for Platelets. *J. Am. Soc. Mass. Spectrom.* 25:651-661
- 18 Sonogo G, Abonnenc M, Tissot JD, Prudent M, Lion N (2017) Redox Proteomics and Platelet Activation: Understanding the Redox Proteome to Improve Platelet Quality for Transfusion. *Int J Mol Sci* 18:387-408
- 19 Lin L, Londe H, Janda JM, Hanson CV, Corash L (1994) Photochemical inactivation of pathogenic bacteria in human platelet concentrates. *Blood* 83:2698-2706

- 20 Lin L, Cook ND, Wieseahn GP et al. (1997) Photochemical inactivation of viruses and bacteria in platelet concentrates by use of a novel psoralen and long-wavelength ultraviolet light. *Transfusion* 37: 423-435
- 21 Lin L (2001) Inactivation of cytomegalovirus in platelet concentrates using Helinx technology. *Semin Hematol* 38:27-33
- 22 Lin L, Dikeman R, Molini B et al. (2004) Photochemical treatment of platelet concentrates with amotosalen and long-wavelength ultraviolet light inactivates a broad spectrum of pathogenic bacteria. *Transfusion* 44:1496-1504
- 23 Lin L, Hanson CV, Alter HJ et al. (2005) Inactivation of viruses in platelet concentrates by photochemical treatment with amotosalen and long-wavelength ultraviolet light. *Transfusion* 45:580-590
- 24 Eastman RT, Barrett LK, Dupuis K, Buckner FS, Van Voorhis WC (2005) Leishmania inactivation in human apheresis platelets by a psoralen (amotosalen HCl) and long-wavelength ultraviolet irradiation. *Transfusion* 45:1459-1463
- 25 Jauvin V, Alfonso RD, Guillemain B, Dupuis K, Fleury HJ (2005) In vitro photochemical inactivation of cell-associated human T-cell leukemia virus Type I and II in human platelet concentrates and plasma by use of amotosalen. *Transfusion* 45:1151-1159
- 26 Pinna D, Sampson-Johannes A, Clementi M, Poli G, Rossini S, Lin L, Vicenzi E (2005) Amotosalen photochemical inactivation of severe acute respiratory syndrome coronavirus in human platelet concentrates. *Transfus Med* 15:269-276
- 27 Gallian P, Vignoli C, Dombey AM et al. (2006) Inactivation of a European strain of West Nile virus in single-donor platelet concentrate using the INTERCEPT™ blood system. *Vox Sang* 91:345-347
- 28 Roback JD, Conlan M, Drew WL, Ljungman P, Nichols WG, Preiksaitis JK (2006) The role of photochemical treatment with amotosalen and UV-A light in the prevention of transfusion-transmitted cytomegalovirus infections. *Transfus Med Rev* 20:45-56
- 29 Brecher ME, Hay S, Corash L, Hsu J, Lin L (2007) Evaluation of bacterial inactivation in prestorage pooled, leukoreduced, whole blood-derived platelet concentrates suspended in plasma prepared with photochemical treatment. *Transfusion* 47:1896-1901
- 30 Nussbaumer W, Allersdorfer D, Grabmer C et al. (2007) Prevention of transfusion of platelet components contaminated with low level bacteria: a comparison of bacteria culture and pathogen inactivation methods. *Transfusion* 47:1125-1133
- 31 Grellier P, Benach J, Labaied M et al. (2008) Photochemical inactivation with amotosalen and long-wavelength ultraviolet light of Plasmodium and Babesia in platelet and plasma components. *Transfusion* 48:1676-1684
- 32 Kwon SY, Kim IS, Bae JE et al. (2014) Pathogen inactivation efficacy of MIRASOL® PRT System and INTERCEPT™ Blood System for non-leucoreduced platelet-rich plasma-derived platelets suspended in plasma. *Vox Sang* 107:254-260
- 33 Musso D, Richard V, Broult J, Cao-Lormeau VM (2014) Inactivation of dengue virus in plasma with amotosalen and ultraviolet A illumination. *Transfusion* 54:2924-2930
- 34 Aubry M, Ricgard V, Green UJ, Broult J, Musso D (2016) Inactivation of Zika virus in plasma with amotosalen und ultraviolet A illumination. *Transfusion* 56:33-40
- 35 CERUS: INTERCEPT™ Plasma Technical Data Sheet (2016)
http://INTERCEPT™bloodsystem.com/sites/default/files/resources/technical_data_sheet_-_plasma_-_english.pdf Zugegriffen 18.05.2018
- 36 CERUS: INTERCEPT™ Platelets Technical Data Sheet (2017)
http://interceptbloodsystem.com/sites/default/files/resources/prdtds_00121_v9.0_secure_0.pdf Zugegriffen 18.05.2018
- 37 Hauser L, Roque-Alphonso AM, Beylouné A et al. (2014) Hepatitis E transmission by transfusion of INTERCEPT™ blood system-treated plasma. *Blood* 123:796-797
- 38 Juhl D, Baylis SA, Blümel J, Görg S, Hennig H (2014) Seroprevalence and incidence of hepatitis E virus infection in German blood donors. *Transfusion* 54:49-56
- 39 Panning M, Umhau M, Emmerich F, Wedemeyer H, Pischke S (2015) Diskussion zu dem Beitrag

- Hepatitis E in Deutschland-eine unterschätzte Infektionskrankheit. *Deutsches Ärzteblatt* 112:220-221
- 40 Baylis SA, Corman VM, Ong E, Linnen JM, Nübling CM, Blümel J (2016) Hepatitis E viral loads in plasma pools for fractionation. *Transfusion* 56:2532-2537
 - 41 Domanovic D, Tedder R, Blümel J et al. (2017) Hepatitis E and blood donation safety in selected European countries: a shift to screening? *Eurosurveillance* 22(16): 30514
 - 42 Schmidt M, Hourfar MK, Sireis W et al. (2015) Evaluation of the effectiveness of a pathogen inactivation technology against clinically relevant transfusion-transmitted bacterial strains. *Transfusion* 55:2104-2112
 - 43 Störmer M, Arroyo A, Brachert J et al. (2012) Establishment of the first international repository for transfusion-relevant bacteria reference strains: ISBT working party transfusion-transmitted infectious diseases (WP-TTID), subgroup on bacteria. *Vox Sang* 102:22-31
 - 44 Benjamin RJ, Wagner SJ (2015) Bacterial pathogen reduction requires validation under conditions of intended use. *Transfusion* 55:2060-2063
 - 45 Grass JA, Hei DJ, Metchette K et al. (1998) Inactivation of leukocytes in platelet concentrates by photochemical treatment with psoralen plus UVA. *Blood* 91:2180-2188
 - 46 Grass JA, Wafa T, Reames A, Wages D, Corsh L, Ferrara JL, Lin L (1999) Prevention of transfusion-associated graft-versus-host disease by photochemical treatment. *Blood* 93:3140-3147
 - 47 Tice RR, Gatehouse D, Kirkland D, Speit G (2007) The pathogen reduction treatment of platelets with S-59 (Amotosalen) plus ultraviolet A light: Genotoxicity profile and hazard assessment. *Mutation Research/Genetic Toxicology and Environmental Mutagenesis* 630:50-68
 - 48 Lin L, Conlan MG, Tessman J, Cimino G, Porter S (2005) Amotosalen interactions with platelet and plasma components: absence of neoantigen formation after photochemical treatment. *Transfusion* 45:1610-1620
 - 49 Van Rhenen DJ, Vermeij J, Mayaudon V, Hind C, Lin L, Corash L (2000) Functional characteristics of S-59 photochemically treated platelet concentrates derived from buffy coats. *Vox Sang* 79:206-214
 - 50 Jansen GAJ, van Vlieth HH, Vermeij H et al. (2004) Functional characteristics of photochemically treated platelets. *Transfusion* 44:313-319
 - 51 Van Aelst B, Feys HB, Devloo R, Vanhoorelbeke K, Vandekerckhove P, Compennolle V (2015) Riboflavin and amotosalen photochemical treatments of platelet concentrates reduce thrombus formation kinetics in vitro. *Vox Sang* 108:328-339
 - 52 Apelseth TO, Bruserud O, Wentzel-Larsen T, Bakken AM, Bjorsvik S, Hervig T (2007) In vitro evaluation of metabolic changes and residual platelet responsiveness in photochemical treated and gamma-irradiated single-donor platelet concentrates during long-term storage. *Transfusion* 47:653-665
 - 53 Johnson L, Loh YS, Kwok M, Marks DS (2013) In vitro assessment of buffy-coat derived platelet components suspended in SSP+ treated with the INTERCEPT™ Blood system. *Transfus Med* 23:121-129
 - 54 Janetzko K, Lin L, Eichler H, Mayaudon V, Flament J, Kluter H (2004) Implementation of the INTERCEPT™ Blood System for Platelets into routine blood bank manufacturing procedures: evaluation of apheresis platelets. *Vox Sang* 86:239-245
 - 55 Picker SM, Speer R, Gathof BS (2004) Functional characteristics of buffy-coat PLTs photochemically treated with amotosalen-HCL for pathogen inactivation. *Transfusion* 44:320-329
 - 56 Sandgren P, Diedrich B (2015) Pathogen inactivation of double-dose buffy-coat platelet concentrates photochemically treated with amotosalen and UVA light: preservation of in vitro function. *Vox Sang* 108:340-349
 - 57 Lozano M, Galan A, Mazzara R, Corash L, Escolar G (2007) Leuco-reduced buffy coat-derived platelet concentrates photochemically treated with amotosalen HCl and ultraviolet A light stored up to 7 days: assessment of hemostatic function under flow conditions. *Transfusion* 47:666-671
 - 58 Picker SM, Schneider V, Gathof BS (2009) Platelet function assessed by shear-induced deposition of split-triple dose apheresis concentrates with pathogen reduction technologies. *Transfusion*

- 49:1224-1232
- 59 Hechler B, Ohlmann P, Chafey P, Ravanat C, Eckly A, Maurer C (2013) Preserved functional and biochemical characteristics of platelet components prepared with amotosalen and ultraviolet A for pathogen inactivation. *Transfusion* 53:1187-1200
 - 60 Tynngard N, Johansson BM, Lindahl TL, Berlin G, Hansson M (2008) Effects of INTERCEPT™ pathogen inactivation on platelet function as analysed by free oscillation rheometry. *Transfus Apher Sci* 38:85-88
 - 61 Osman A, Hitzler WE, Meyer CU et al. (2015) Effects of pathogen reduction systems on platelet microRNAs, mRNAs, activation, and function. *Platelets* 26:154–163
 - 62 Bakkour S, Chafets DM, Wen L et al. (2016) Assessment of nucleic acid modification induced by amotosalen and ultraviolet A light treatment of platelets and plasma using real-time polymerase chain reaction amplification of variable length fragments of mitochondrial DNA. *Transfusion* 56:410-420.
 - 63 Laffont B, Corduan A, Ple H et al. (2013) Activated platelets can deliver mRNA regulatory Ago2 microRNA complexes to endothelial cells via microparticles. *Blood* 122:253-261
 - 64 Semple JW (2013) Platelets deliver small packages of genetic function. *Blood* 122:155-156
 - 65 Osman A, Hitzler WE, Ameer A, Provost P (2015) Differential Expression Analysis by RNA-Seq Reveals Perturbations in the Platelet mRNA Transcriptome Triggered by Pathogen Reduction Systems. *PlosOne* DOI:10.1371/journal.pone.0133070
 - 66 Prudent M, Crettaz D, Delobel J, Tissot JD, Lion N J (2012) Proteomic analysis of Intercept-treated platelets. *Proteomics* 76:316-328
 - 67 Prudent M, D'Alessandro A, Cazenave JP et al. (2014) Proteome changes in platelets after pathogen inactivation--an interlaboratory consensus. *Transfus Med Rev.* 28:72-83
 - 68 Terumo BCT MIRASOL® Pathogen Reduction Technology System, 05.12.2014. <https://www.terumobct.com/location/north-america/products-and-services/Pages/MIRASOL®-pathogen-reduction-technology.aspx> Zugriffen: 22.05.2018
 - 69 Corbin F: Pathogen inactivation of blood components (2002) current status and introduction of an approach using riboflavin as a photosensitizer. *Int J Hematol* 76 Suppl 2:253-257
 - 70 Marschner S, Goodrich R (2011) Pathogen reduction technology treatment of platelets, plasma and whole blood using riboflavin and UV light. *Transfus Med Hemother* 38:8-18
 - 71 Owada T, Kaneko M, Matsumoto C et al. (2014) Establishment of culture systems for Genotypes 3 and 4 hepatitis E virus (HEV) obtained from human blood and application of HEV inactivation using a pathogen reduction technology system. *Transfusion* 54(11):2820-7
 - 72 Ruane PH, Edrich R, Gampp D, Keil SD, Leonard RL, Goodrich RP (2004) Photochemical inactivation of selected viruses and bacteria in platelet concentrates using riboflavin and light. *Transfusion* 44:877-885
 - 73 Keil SD, Bengrine A, Bowen R et al. (2015) Inactivation of viruses in platelet and plasma products using a riboflavin-and-UV-based photochemical treatment. *Transfusion* 55:1736-1744
 - 74 Faddy HM, Fryk JJ, Watterson D et al. (2016) Riboflavin and ultraviolet light impact on dengue virus infectivity. *Vox Sang* 111:235-241
 - 75 Faddy HM, Prow NA, Fryk JJ et al. (2015) The effect of riboflavin and ultraviolet light on the infectivity of arboviruses. *Transfusion* 55:825-831
 - 76 Vanlandingham DL, Keil SD, Horne KM, Pyles R, Goodrich RP, Higgs S (2013) Photochemical inactivation of chikungunya virus in plasma and platelets using the MIRASOL® pathogen reduction technology system. *Transfusion* 53:284-290
 - 77 Keil SD, Saakadze N, Bowen R et al. (2015) Riboflavin and ultraviolet light for pathogen reduction of murine cytomegalovirus in blood products. *Transfusion* 55:858-863
 - 78 Goodrich RP, Edrich R, Li J, Segatchian J (2006) The MIRASOL® PRT system for pathogen reduction of platelets and plasma: an overview of current status and future trends. *Transfus Apher Sci* 35:5-17
 - 79 Shinohara N, Matsumoto C, Chatani M et al. (2015) Efficacy of the MIRASOL® pathogen reduction technology system against severe fever with thrombocytopenia syndrome virus (SFTSV). *Vox*

- Sang 109:417-419
- 80 Rentas F, Harman R, Gomez C et al. (2007) Inactivation of *Orientia tsutsugamushi* in red blood cells, plasma, and platelets with riboflavin and light, as demonstrated in an animal model. *Transfusion* 47:240-247
 - 81 Tonnetti L, Proctor MC, Reddy HL, Goodrich RP, Leiby DA (2010) Evaluation of the MIRASOL® pathogen reduction technology system against *Babesia microti* in apheresis platelets and plasma. *Transfusion* 50:1019-1027
 - 82 Tonnetti L, Thorp AM, Reddy HL, Keil SD, Goodrich RP, Leiby DA (2013) Riboflavin and ultraviolet light reduce the infectivity of *Babesia microti* in whole blood. *Transfusion* 53:860-867
 - 83 Tonnetti L, Thorp AM, Reddy HL, Keil SD, Goodrich RP, Leiby DA (2015) Reduction of *Leishmania donovani* infectivity in whole blood using riboflavin and ultraviolet light. *Transfusion* 55:326-329
 - 84 Chaar ME, Atwal S, Freimanis GL, Dinko B, Sutherland CJ, Allain J-P (2013) Inactivation of *Plasmodium falciparum* in whole blood by riboflavin plus irradiation. *Transfusion* 53:3174-3183
 - 85 Keil SD, Kiser P, Sullivan SS et al. (2013) Inactivation of plasmodium spp. In plasma and platelet concentrates using riboflavin and ultraviolet light. *Transfusion* 53:2278-2286
 - 86 Tonnetti L, Thorp AM, Reddy HL, Keil SD, Goodrich RP, Leiby DA (2012) Evaluating pathogen reduction of *Trypanosoma cruzi* with riboflavin and ultraviolet light for whole blood. *Transfusion* 52:409-416
 - 87 Cardo LJ, Salata J, Mendez J, Reddy H, Goodrich R (2007) Pathogen inactivation of *Trypanosoma cruzi* in plasma and platelet concentrates using riboflavin and ultraviolet light. *Transfus Apher Sci* 37:131-137
 - 88 Brecher ME, Hay SN (2005) Bacterial contamination of blood components. *Clin Microbiol Rev* 18: 195-204
 - 89 Goodrich RP, Gilmour D, Hovenga N, Keil SD (2009) A laboratory comparison of pathogen reduction technology treatment and culture of platelet products for addressing bacterial contamination concerns. *Transfusion* 49:1205-1216
 - 90 Fast LD, DiLeone G, Li J, Goodrich RP (2006) Functional inactivation of white blood cells by MIRASOL® treatment. *Transfusion* 46:642-648
 - 91 Fast LD, DiLeone G, Cardarelli G, Li J, Goodrich RP (2006) MIRASOL® PRT treatment of donor white blood cells prevents the development of xenogeneic graft-versus-host disease in Rag2^{-/-}γc^{-/-} double knockout mice. *Transfusion* 46:1553-1560
 - 92 Jackman RP, Heitman JW, Marschner S, Goodrich RP, Norris PJ (2009) Understanding loss of donor white blood cell immunogenicity after pathogen reduction: mechanisms of action in ultraviolet illumination and riboflavin treatment. *Transfusion* 49:2686-2699
 - 93 Marschner S, Fast LD, Baldwin III WM, Slichter SJ, Goodrich RP (2010) White blood cell inactivation after treatment with riboflavin and ultraviolet light. *Transfusion* 50:2489-2498
 - 94 Fast LD, DiLeone G, Marschner S (2011) Inactivation of human white blood cells in platelet products after pathogen reduction technology treatment in comparison to gamma irradiation. *Transfusion* 51:1397-404
 - 95 Fast LD, Nevala M, Tavares J, Reddy H, Goodrich RP, Marschner S (2013) Treatment of whole blood with riboflavin plus ultraviolet light, an alternative to gamma irradiation in the prevention of transfusion-associated graft-versus-host disease? *Transfusion* 53:373-381
 - 96 Reddy HL, Doane SK, Keil SK, Marschner S, Goodrich RP (2013) Development of a riboflavin and ultraviolet light-based device to treat whole blood. *Transfusion* 53:131S-136S
 - 97 FDA: Database of Select Committee on GRAS Substances (SCOGS) Reviews: Riboflavin <https://www.fda.gov/Food/IngredientsPackagingLabeling/GRAS/SCOGS/default.htm> Zugegriffen 22.05.2018
 - 98 Diakonou VF, Grentzelos MA, Tzatzarakis MN et al. (2012) Riboflavin's time-dependent degradation rate induced by ultraviolet A irradiation. *Eur J Ophthalmol* 22 Suppl. 7: 51-56
 - 99 Ambruso DR, Thurman G, Marschner S, and Raymond P. Goodrich RP (2009) Lack of antibody formation to platelet neoantigens after transfusion of riboflavin and ultraviolet light-treated platelet concentrates. *Transfusion* 49:2631-2636

- 100 Janetzko K, Hinz K, Marschner S, Goodrich R, Klüter H (2009) Pathogen reduction technology (MIRASOL®) treated single-donor platelets resuspended in a mixture of autologous plasma and PAS. *Vox Sang* 97:234-239
- 101 Li J, Goodrich L, Hansen E, Edrich R, Gampp D, Goodrich RP (2005) Platelet glycolytic flux increases stimulated by ultraviolet-induced stress is not the direct cause of platelet morphology and activation changes: possible implications for the role of glucose in platelet storage. *Transfusion* 45:1750-1758
- 102 Picker SM, Tauszig ME, Gathof BS (2012) Cell quality of apheresis-derived platelets treated with riboflavin-ultraviolet light after resuspension in platelet additive solution. *Transfusion* 52:510-516
- 103 Ostrowski SR, Bochsén L, Salado-Jimena JA et al. (2010) In vitro cell quality of buffy coat platelets in additive solution treated with pathogen reduction technology. *Transfusion* 50:2210-2219.
- 104 Johnson L, Winter KM, Reid S et al. (2011) The effect of pathogen reduction technology (MIRASOL®) on platelet quality when treated in additive solution with low plasma carryover. *Vox Sang* 101:208–214
- 105 Cookson P, Thomas S, Marschner S, Goodrich R, Cardigan R (2012) In vitro quality of single-donor platelets treated with riboflavin and ultraviolet light and stored in platelet storage medium for up to 8 days. *Transfusion* 52:983-994
- 106 Perez-Pujol S, Tonda R, Lozano M et al. (2005) Effects of a new pathogen-reduction technology (MIRASOL® PRT) on functional aspects of platelet concentrates. *Transfusion* 45:911-919
- 107 Picker SM, Steisel A, Gathof BS (2008) Effects of MIRASOL® PRT treatment on storage lesion development in plasma-stored apheresis-derived platelets compared to untreated and irradiated units. *Transfusion* 48:1685-1692
- 108 Picker SM, Steisel A, Gathof BS (2009) Evaluation of White Blood Cell- and Platelet-Derived Cytokine Accumulation in MIRASOL®-PRT-Treated Platelets. *Transfus Med Hemother* 36:114–120
- 109 Picker SM, Oustianskaia L, Schneider V, Gathof BS (2010) Annexin V Release and Transmembrane Mitochondrial Potential during Storage of Apheresis-Derived Platelets Treated for Pathogen Reduction. *Transfus Med Hemother* 37:7–12
- 110 Castrillo A, Cardoso M, Rouse L (2013) Treatment of Buffy Coat Platelets in Platelet Additive Solution with the MIRASOL® Pathogen Reduction Technology System. *Trans Med Hemother* 40:44–48
- 111 Zeddies S, De Cuyper IM, van der Meer PF et al. (2014) Pathogen reduction treatment using riboflavin and ultraviolet light impairs platelet reactivity toward specific agonists in vitro. *Transfusion* 54:2292-2300
- 112 Galan AM, Lozano M, Molina P et al. (2011) Impact of pathogen reduction technology and storage in platelet additive solutions on platelet function. *Transfusion* 51:808-815
- 113 Terada C, Mori J, Okazaki H, Satake M, Tadokoro K (2014) Effects of riboflavin and ultraviolet light treatment on platelet thrombus formation on collagen via integrin $\alpha\text{IIb}\beta\text{3}$ activation. *Transfusion* 54:1808-1818
- 114 Klein-Bosgoed C, Schubert P, Devine DV (2016) Riboflavin and ultraviolet illumination affects selected platelet mRNA transcript amounts differently. *Transfusion* 56:2286-2295
- 115 Feys HB, van Aelst B, Devreese K et al. (2014) Oxygen removal during pathogen inactivation with riboflavin and UV preserves protein function in plasma for transfusion. *Vox Sang* 106:307-315
- 116 Johnson L, Marks D (2015) Treatment of Platelet Concentrates with the MIRASOL® Pathogen Inactivation System Modulates Platelet Oxidative Stress and NF- κB Activation. *Transfus Med Hemother* 42:167–173
- 117 Schubert P, Culibrk B, Coupland D, Scammell K, Gyongyossy-Issa M, Devine DV (2012) Riboflavin and ultraviolet light treatment potentiates vasodilator-stimulated phosphoprotein Ser-239 phosphorylation in platelet concentrates during storage. *Transfusion* 52:397–408
- 118 AuBuchon JP, Herschel L, Roger J et al. (2005) Efficacy of apheresis platelets treated with riboflavin and ultraviolet light for pathogen reduction. *Transfusion* 45:1335–1341.
- 119 Goodrich RP, Li J, Pieters H, Crookes R, Roodt J, Heyns AdP (2006) Correlation of in-vitro platelet quality measurements with in vivo platelet viability in human subjects. *Vox Sang* 90:279–285
- 120 Macopharma Theraflex UV-Platelets (2014) <http://blood->

- safety.macopharma.com/en/category/documents-literatures/theraflex-uv-platelets Zugegriffen 22.05.2018
- 121 Seltsam A, Muller TH (2011) UVC Irradiation for Pathogen Reduction of Platelet Concentrates and Plasma. *Transfus Med Hemother* 38:43-54.
 - 1222 Seghatchian J, Tolksdorf F (2012) Characteristics of the THERAFLEX UV-Platelets pathogen inactivation system - an update. *Transfus Apher Sci* 46:221-229
 - 123 Mohr H, Steil L, Gravemann U et al. (2009) A novel approach to pathogen inactivation in platelet concentrates using short wave ultraviolet light. *Transfusion* 49:2612-2624
 - 124 Thiele T, Pohler P, Kohlmann T et al. (2015) Tolerance of platelet concentrates treated with UVC-light only for pathogen reduction – a phase I clinical trial. *Vox Sang* 109:44-51
 - 125 Kim S, Handke W, Gravemann U, Döscher A, Brixner V, Müller TH, Seltsam A (2018) Mitochondrial DNA multiplex real-time polymerase chain reaction inhibition assay for quality control of pathogen inactivation by ultraviolet C light in platelet concentrates *Transfusion* DOI: 10.1111/trf.14464
 - 126 Steinmann E, Gravemann U, Friesland M, Dörrbecker J, Müller T H, Pietschmann T, Seltsam A (2013) Two pathogen reduction technologies--methylene blue plus light and shortwave ultraviolet light--effectively inactivate hepatitis C virus in blood products. *Transfusion* 53: 1010-1018
 - 127 van der Meer P F, Gravemann U, de Korte D, Sumian C, Tolksdorf F, Muller TH, Seltsam A (2016) Effect of increased agitation speed on pathogen inactivation efficacy and in vitro quality in UVC-treated platelet concentrates. *Vox Sang* 111: 127-134
 - 128 Faddy HM, Fryk JJ, Prow NA et al. (2016) Inactivation of dengue, chikungunya, and Ross River viruses in platelet concentrates after treatment with ultraviolet C light. *Transfusion* 56: 1548-1555
 - 129 Seltsam A (2011) Pathogenreduktion mittels UVC-Licht: Entwicklung des THERAFLEX UV-Plättchen-Systems. DGTI, Freiburg
<https://www.dgti.de/docs/doclink/10802/6.%20Vortrag%20Seltsam,%20Theraflex%20UV%20Platelets.pdf>
 Zugegriffen 23.05.2018
 - 130 Gravemann U, Schmidt JP, Tolksdorf F, Sumian C, Müller TH, Bauer G, Biczysko S, Keppler OT, Seltsam A, Rabenau HF (2013) Influenza A virus H3N2 is efficiently inactivated by the THERAFLEX UV-PLATELETS system. *Vox Sang* 105 (Suppl. 1), 134
 - 131 Fryk JJ, Marks DC, Hobson-Peters J, Watterson D, Hall RA, Young PR, Reichenberg S, Tolksdorf F, Sumian C, Gravemann U, Seltsam A, Faddy HM (2017) Reduction of Zika virus infectivity in platelet concentrates after treatment with ultraviolet C light and in plasma after treatment with methylene blue and visible light. *Transfusion* 57:2677-2682
 - 132 Gravemann U, Lambrecht B, Schmidt J-P, Seltsam A (2014) Hepatitis A virus is efficiently inactivated by the THERAFLEX UV-Platelets system. *Vox Sang* 107 (Suppl. 1), 128
 - 133 Spindler-Raffel E, Benjamin RJ, McDonald CP, Ramirez-Arcos S, Aplin K, Bekeredjian-Ding I, de Korte D, Gabriel C, Gathof B, Hanschmann KM, Hourfar K, Ingram C, Jacobs MR, Keil SD, Kou Y, Lambrecht B, Marcelis J, Mukhtar Z, Nagumo H, Niekerk T, Rojo J, Marschner S, Satake M, Seltsam A, Seifried E, Sharafat S, Stormer M, Sussner S, Wagner SJ, Yomtovian R. ISBT Working Party Transfusion-Transmitted Infectious Diseases, Subgroup on Bacteria (2017) Enlargement of the WHO international repository for platelet transfusion-relevant bacteria reference strains. *Vox Sang* 112: 713-722
 - 134 Gravemann U, Tolksdorf F, Handke W, Müller TH, Seltsam A (2017) The THERAFLEX UV-Platelets technology efficiently inactivates transfusion-relevant bacteria species in contaminated platelet concentrates. *Transfusion* 57 (Suppl. 1), 203A
 - 135 Castro, E, Girones N, Guerrero N, Barea L, Fresno M (2008) The effectiveness of UVC pathogen inactivation system on reducing the *Trypanosoma cruzi* and *Leishmania infantum* burden in platelets. *Vox Sang* 95 (Suppl. 1), 290
 - 136 Castro E, González LM, Bautista JM, Montero E, Rubio JM (2013) The efficacy of UVC pathogen inactivation on the reduction of *Plasmodium falciparum* in buffy coat derived platelets. *Vox Sang*

- 105 (Suppl. 1), 133
- 137 Castro E, González LM, Rubio JM, Ramiro R, Gironés N, Montero E (2014) The efficacy of the ultraviolet C pathogen inactivation system in the reduction of *Babesia divergens* in pooled buffy coat platelets. *Transfusion* 54: 2207-2216
- 138 Pohler P, Müller M, Winkler C et al. (2015) Pathogen reduction by ultraviolet C light effectively inactivates human white blood cells in platelet products. *Transfusion* 55: 337-347
- 139 Pohler P, Lehmann J, Veneruso V et al. (2012) Evaluation of the tolerability and immunogenicity of ultraviolet C-irradiated autologous platelets in a dog model. *Transfusion* 52:2414-2426
- 140 Bashir S, Cookson P, Wiltshire M et al. (2013) Pathogen inactivation of platelets using ultraviolet C light effect on in vitro function and recovery and survival of platelets. *Transfusion* 53:990-1000
- 141 Tynngard N, Trinks M, Berlin (2015) In vitro function of platelets treated with ultraviolet C light for pathogen inactivation: a comparative study with non-irradiated and gamma irradiated platelets. *Transfusion* 55:1169–1177
- 142 Van Aelst B, Devloo R, Vandekerckhove P, Compernelle V, Feys HB (2015) Ultraviolet C light pathogen inactivation treatment of platelet concentrates preserves integrin activation but affects thrombus formation kinetics on collagen in vitro. *Transfusion* 55:2404-2414
- 143 Martini R, Kempfer CB, Rodrigues Mde A et al. (2010) Bacterial contamination on platelet concentrates: identification, antimicrobial susceptibility profile and transfusion-related sepsis. *Revista da Sociedade Brasileira de Medicina Tropical* 43:682-685
- 144 Kou Y, Pagotto F, Hannach B, Ramirez-Arcos S (2015) Fatal false-negative transfusion infection involving a buffy coat platelet pool contaminated with biofilm-positive *Staphylococcus epidermidis*: a case report. *Transfusion* 55:2384-2389
- 145 Snyder E, Raife T, Lin L et al. (2004) Recovery and life span of 111indium-radiolabeled platelets treated with pathogen inactivation with amotosalen HCl (S-59) and ultraviolet A light. *Transfusion* 44:1732–1740
- 146 Murphy S (2004) Radiolabeling of PLTs to assess viability: a proposal for a standard. *Transfusion* 44:131-133.
- 147 FDA Blood Products Advisory Committee (2004) July 22-23th, Summary Topic II: New Standards for Platelet Evaluation. <https://wayback.archive-it.org/7993/20170404063056/https://www.fda.gov/ohrms/dockets/ac/04/transcripts/2004-4057t1.htm> Zugegriffen 23.05.2018
- 148 Swiss Medic: Haemovigilance Annual report 2015 www.swissmedic.ch/haemovigilance-report Zugegriffen 23.05.2018
- 149 Butler C, Doree C, Estcourt LJ et al. (2013) Pathogen-reduced platelets for the prevention of bleeding. *Cochrane Database Syst Rev.* 28:CD009072
- 150 Estcourt LJ, Malouf R, Hopewell S et al. (2017) Pathogen-reduced platelets for the prevention of bleeding. *Cochrane Database Syst Rev.* Issue 7. Art. No.: CD009072
- 151 McCullough J, Vesole DH, Benjamin RJ et al. (2004) Therapeutic efficacy and safety of platelets treated with a photochemical process for pathogen inactivation: the SPRINT Trial. *Blood* 104:1534-1541
- 152 Kerkhoffs JLH, van Putten WLJ, Novotny VMJ et al. (2010) Clinical effectiveness of leucoreduced, pooled donor platelet concentrates, stored in plasma or additive solution with and without pathogen reduction. *BJH* 150:209–217
- 153 Cid J, Escolar G, Lozano M (2012) Therapeutic efficacy of platelet components treated with amotosalen and ultraviolet A pathogen inactivation method: results of a meta-analysis of randomized controlled trials. *Vox Sang* 103:322-330
- 154 Garban F, Guyard A, Labussière H et al. (2018) Comparison of the Hemostatic Efficacy of Pathogen-Reduced Platelets vs Untreated Platelets in Patients With Thrombocytopenia and Malignant Hematologic Diseases: A Randomized Clinical Trial. *JAMA Oncol* 4(4):468-475
- 155 The MIRASOL® Clinical Evaluation Study Group (2010) A randomized controlled clinical trial evaluating the performance and safety of platelets treated with MIRASOL®. *Transfusion* 50: 2362-2375

- 1565 Snyder E, McCullough J, Slichter SJ et al. (2005) Clinical safety of platelets photochemically treated with amotosalen HCl and ultraviolet A light for pathogen inactivation: the SPRINT trial. *Transfusion* 45: 1864-1875
- 1576 Kleinman S, Reed W, Stassinopoulos A (2013) A patient-oriented risk-benefit analysis of pathogen-inactivated blood components: application to apheresis platelets in the United States. *Transfusion* 53: 1603-1618
- 158 Amato M, Schennach H, Astl M et al. (2017) Impact of platelet pathogen inactivation on blood component utilization and patient safety in a large Austrian Regional Medical Centre. *Vox Sang* 112:47-55
- 159 Knutson F, Osselaer J, Pierelli L et al. (2015) A prospective, active haemovigilance study with combined cohort analysis of 19 175 transfusions of platelet components prepared with amotosalen-UVA photochemical treatment. *Vox Sang* 109:343-52
- 160 Jutzi M, Rüesch M, Taleghani BM (2013) Einführung der Pathogeninaktivierung für Thrombozytenkonzentrate in der Schweiz. *Schweiz Med Forum* 13:222–226

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